ANIMAL & CELL HYPOXIC WORKSTATION

This unique workstation can be used for hypoxia cells and animal research. Two workstations are flexible and divided, and each workstation is controlled independently.



Double workstation: the left workstation is used for the hypoxia study of animals, and the right workstation is used for the hypoxia study of cells.

Operation	Cultivation, operation and observation in one place Provides stable temperature, humidity, and gas control within the workstation			
	• Easy to operate Direct access to the workstation in only 8-10 seconds via vacuum Avoid cross-contamination and skin allergies caused by traditional rubber gloves			
Oystem	 Accurate gas control Real-time monitoring and control the O₂, CO₂ concentration, temperature and humidity 			
	• Efficiency & Reliability The standard built-in interlock allows easy transfer of animals or supplies into the workstation			
	 Large Operating Space The DY-WC-D workstation has two spacious areas that allow a variety of monitoring or imaging systems to be integrated inside the workstation, as well as a large number of samples 			
	 • O₂ control range: 0.5% - 20.9%, in increments of 0.1% 			
	 CO₂ alarm system: Avoid high-concentration carbon dioxide that affects the experiment 			
Specifications (Left)	 Programmable four-stage hypoxia cycle: Allows a user-defined timed sequence of up to 4 different O₂ concentrations 			
	 Temperature range: 5°C above ambient up to 45°C, in 0.1°C increments 			
	 Data recording: USB data logging system with continuous data storage 			

· 12 levels of adjustable light intensity to simulate normal sunlight exposure

	• O ₂ control range: 0.1% - 23.0%, in increments of 0.1%			
	CO ₂ control range: 0.03% - 20.0%, in increments of 0.1%			
Specifications	• Programmable four-stage hypoxia cycle: Allows a user-defined timed sequence of up to 4 different O ₂ and CO ₂ concentrations			
(Right)	• Temperature range: 5°C above ambient up to 45°C , in 0.1°C increments			
	Humidity range: Room humidity to 75%, in 0.1% increments			
	 One-touch calibration for O₂ and CO₂ sensor 			
• Data recording: USB data logging system with 3 months of continuous data storage				
Standard System	 Direct-access operation system Cylinder low pressure alarm USB data logging system Internal socket Sensor calibration Interlock purge at different O₂ Levels 	 Internal light Activated carbon filtration system Humidity control (R) Temperature control (R) Desiccants (L) 		
Optional System	 Foot Pedal Interlock Purge Four-stage hypoxia cycle HEPA system Removable front panel Cable gland Vacuum Suction Pipe 	 Single plate entry system (R) Anaerobic Mode (R) Temperature control (L) CO₂ Adsorbent (L) CO₂ alarm system (L) 12 levels internal light (L) 		

Application Field sleep apnoea syndrome, hyperoxia hypoxia, pulmonary hypertensionNeurological disorders: cerebral ischaemic injury, neurobiology

· Cardiovascular: myocardial ischaemia and hypoxia

Oncology: tumour biology

· Respiratory related:

Model		L	R
Internal size	Height (mm)	420	550
	Width(mm)	800	760
	Depth(mm)	460	570
Interlock size	Height (mm)	245	290
	Width(mm)	288	270
	Depth(mm)	257	345

MAWORDE

Tel: 8610-88693537 Email: sales@maworde.com www.maworde-biotech.com