

ANIMAL INFECTION INCUBATOR

According to the set gas concentration, the proportioning gas can be connected to a specific gas, such as smog/toxic gas for pathological research. Or connect O₂, N₂, CO₂ for related hypoxia research.



Control System

- Automatic control system mixed gas enters the experiment box and controls the internal pressure of the box
- Gas control range: 0.5% - 20.9%, in increments of 0.1%
- Optional temperature control: 5°C above ambient up to 45°C , in 0.1°C increments
- Can achieve low oxygen or a specific gas environment inside the chamber
- Customize the gas intake time and the interior of the cabinet

Application Field

- Respiratory related: sleep apnoea syndrome, hyperoxia hypoxia, pulmonary hypertension
- Neurological disorders: cerebral ischaemic injury, neurobiology
- Cardiovascular: myocardial ischaemia and hypoxia
- Oncology: tumour biology

Tel: 8610-88693537

Email: sales@maworde.com

www.maworde-biotech.com

